Ventilateur **OSIRIS**

YL024400 / Rev.2.1 / 10-2012

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FR

User manual

EN

Manual de uso

ES

Manual d'uso

IT

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PL

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PT

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GENERAL SAFETY INSTRUCTIONS

USING OXYGEN

Safety instructions relative to oxygen must be respected:

- No smoking,
- No naked flames or incandescent objects,
- No greasy materials to be used.

Electrical supply

- Check that the mains supply to be used corresponds exactly to the electrical characteristics of the equipment.
- Before the initial use of the equipment and/or when replacing the battery-pack, it is necessary to completely charge the battery-pack (about 10 hours).
- This pulmonary ventilator must not be used with inflammable anaesthetic agents.
- This equipment conforms to the agreed standard EN 60601-1-2 relative to electromagnetic compatibility of medical equipment and to special limitations of EN 794-3 relative to emergency ventilators, but may be affected by adjacent use of apparatus for diathermy, high-frequency electro-surgery, defibrillators or short wave therapy equipment, cellular telephones, or general electromagnetic interference above levels fixed by the EN 60 601-1-2 and EN 794-3 standards.
- Do not use **OSIRIS 2** in specifically magnetic environments (IRM, etc.).
- This equipment conforms to protection requirements of the 93/42/CEE directive.
- A new battery will only provide maximum performance when it has been charged and discharged two or three times!

Gas supplies to the equipment

- Supply the equipment from wall outlets or bottles at a pressure in the range 280 kPa (2.8 bar) to 600 kPa (6 bar).

Starting up the equipment

- Each time the equipment is used check that the audible and visible alarms are working correctly, and go through the check list in the Appendix (page 19).

Use and maintenance of the equipment

Conforming to the standard EN 60 601-1 (Appendix A § 6.8.2.b):

«The manufacturer, the assembler, the installer, or the importer can only be considered responsible for safety aspects, reliability, or the characteristics of the equipment if:

- the assembly, additions, adjustments, modifications, or repairs have been made by persons authorized by him,
- the electrical installation where the equipment is used conforms to CEI instructions,
- the equipment is used according to its instructions for use».

If spare parts used by a technician for periodical maintenance do not conform to manufacturing requirements, the manufacturer is relieved of all responsibility in the case of any incident.

- Do not open a unit which is connected to the mains supply.
- This pulmonary ventilator must not be used with inflammable anaesthetic agents or explosive products..
- Do not use piping or tubing that is antistatic or that conducts electricity.
- **OSIRIS 2** must not be used on a patient without surveillance.
- Do not use **OSIRIS 2** in a hyperbaric chamber.
- For **OSIRIS 2** to function correctly, keep the rear air intakes free of any obstacles, to allow free air circulation.

- In the event of operation in dangerous or explosive atmospheres, it is imperative that **OSIRIS 2** is switched to position VAC/ACMV 100 % O₂.
- The manufacturer has foreseen most cases of possible malfunction and these are normally covered by the integrated alarm system. However, it is recommended that for totally dependent patients, back-up equipment is closely available, such as a manual resuscitator (IM5, etc.).
- Efficient ventilation can be ensured by good settings of the threshold for patients' breath detection (for restless patients, etc.) as well as good surveillance of the patient circuit.
- Using other equipment which does not conform to the safety standards of **OSIRIS 2** can cause a reduction in safety levels. The factors to take into account when choosing other equipment are its use close to the patient, conformity to the EN 794-3 standard and recommendations contained in the accompanying documentation.

Non-invasive ventilation

- Perfect knowledge of **OSIRIS 2** and of respiratory physiopathology is necessary to ensure efficient mask ventilation for the benefit of the patient.



Because of leaks (sometimes unavoidable) around a mask, inhaled volumes may be appreciably different from set volumes.

For certain types of use, it is recommended that **OSIRIS 2** be used with currently approved monitoring equipment.

Equipment cleaning and maintenance

- Do not use abrasive powders, alcohol, acetone, or other easily inflammable solvents.
- **OSIRIS 2** must be checked regularly. To plan and record maintenance operations, refer to the maintenance form in the Appendix.
- Only people who have completely read and understood this manual are authorized to handle and use **OSIRIS 2**.
- The present manual was written to give all the information necessary for using **OSIRIS 2**, but is no substitute for medical instructions, indispensable for adjustments adapted to the needs of the patient.
- The audible alarm of **OSIRIS 2** is intended to be heard by staff close to the patient. Users must determine maximum safe distance from the patient to be able to hear the alarm clearly in the ambient sound environment.

GENERAL SAFETY INSTRUCTIONS	1
I. INTRODUCTION	3
II. DESCRIPTION AND SETTINGS	3
II. 1. Front panel	3
II. 2. Electrical and pneumatic connections on side of unit	5
II. 3. Rear of unit.....	6
III. SETTING UP THE UNIT	7
III. 1. Installing the battery	7
III. 2. Connections.....	7
III. 3. Adjusting OSIRIS 2.....	9
<i>III. 4. DISPLAY OF FREQUENCY AND I/E RATIO MEASURED BY OSIRIS 2.....</i>	<i>12</i>
IV. DO NOT FORGET	13
IV. 1. Safety and protection.....	13
IV. 2. Alarms	13
IV. 3. Operating incidents.....	14
IV. 4. Cleaning, Disinfection, Sterilization	15
IV. 5. Maintenance.....	16
IV. 6. Accessories	17
IV. 7. Waste disposal	17
IV. 8. Disposal of equipment.....	17
V. TECHNICAL DESCRIPTION	18
V. 1. Operational principles	18
V. 2. Technical characteristics.....	18
V. 3. STANDARDS	18
V. 4. Symbols and abbreviations.....	19
VI. APPENDIX.....	20
Material in direct or indirect contact with the patient	20
Pneumatic operation schematic	20
CHECK-LIST: Before use / Return from maintenance.....	21
Maintenance form.....	22

I. INTRODUCTION

I. 1. Description of the Device

OSIRIS 2 is an adjustable, controlled frequency ventilator, with an adjustable I/E ratio from 1/1 to 1/3 by steps of 0.5.

It ensures controlled volume ventilation with either:

- pure oxygen (VAC/ACMV 100 % O₂),
- or with a mixture of O₂ / Ambient air (FiO₂ = about 60 %) (VAC/ACMV Air + O₂).

It also ensures spontaneous ventilation, pressure assisted (VS-AI/PS) which may be made non-invasively.

The PEEP* function is integrated and adjustable from 0 to 15 mbar (0 to 15 cmH₂O).

OSIRIS 2 has audible and visual alarms to signal any ventilation faults. The setting of a maximum pressure threshold (Pmax) fixes a pressure limit and therefore avoids any barotrauma (insufflation stops as soon as the set threshold is passed, and the audible and visual alarms are activated).

PEEP* : Positive End Expiratory Pressure

I. 2. Intended Use

OSIRIS 2 is intended for use by hospital personnel (physicians, nurses, etc.) and non-hospital emergency personnel. It is used:

- For transportation between hospitals,
- For emergency transport between hospitals,
- In post-operative recovery rooms,
- In intensive care units,
- For resuscitation.

For post-operative or intensive care use, **OSIRIS 2** must have currently approved monitoring equipment.

I. 3. Composition of OSIRIS 2

- The basic **OSIRIS 2** unit .
- Three complete single-use patient circuits.
- A 3 meter O₂ supply tube.
- One adult mask.
- One child's mask.
- A mains supply.
- A carrying case.

II. DESCRIPTION AND SETTINGS

II. 1. FRONT PANEL

Identification of parts

① Display screen - parameters, settings, measurements. 2-line display.

② Ventilation mode switch


This switch allows the unit to be made operational. It must be switched to one of the 3 positions: **VAC/ACMV 100%** or **VAC/ACMV Air + O₂** or **VS-AI/PS**.

To stop the unit, turn the switch to position 0.

③ PEEP setting

This knob allows the PEEP to be set, from 0 to 15 cmH₂O.

A safety device ②⑤ must be released to give access to the extended range of 10 to 15 cmH₂O.


④ Safety pressure setting (Pmax)

To make the adjustment:

- Close off the patient output.
- Read the maximum pressure reached in the patient circuit (on the bargraph at the end and after inhalation).
- Adjust Pmaxi alarm threshold by turning knob ④ and check concordance (to within 5 cm H₂O) between the value set with knob ④ (graduated) and the bargraph value.

⑤ Inhibit audible alarm button

The audible alarm can be inhibited for between 1 and 2 minutes by pressing on the  button, but the red lamp will continue to flash and a pilot lamp (5b) will indicate the audible alarm inhibition.

⑥ Battery lamp

This green lamp **flashes** when **OSIRIS 2** is operating **normally**, running on its internal battery pack.

- Short pulse flashing: indicates that OSIRIS 2 is working on batteries.
- Steady flashing: indicates that the battery is charging slowly.
- Constant light: indicates that the battery is charging rapidly.
- Indicator off: Charging complete (charge care).

⑦ Weak battery pack indicator (orange)

This lights up when the level of charge of the internal battery pack has fallen to a point where autonomy of the **OSIRIS 2** is less than 60 minutes..

When this lamp lights up, it is advisable to connect directly to the mains supply or to change the battery pack. Battery recharging is necessary. To charge the battery, plug into the mains supply for 10 hours, with **OSIRIS 2** switched off.

⑧ Mains supply lamp (green)

This is lit when **OSIRIS 2** is connected to a mains supply.

⑨ Frequency setting knob

This knob is used to set minimum frequency cycles to between 6 and 40 c/min.

⑩ Knob for setting inspiratory volume or pressure support


This knob is for setting flowrate in VAC/ACMV 100% O₂ mode (or in VAC/ACMV Air + O₂ mode). The set value is displayed on line 2 of the screen ①.


⑪ "☰" key

This key allows several parameters to be selected for adjustment: I/E ratio, Inspiratory trigger, Pmini and Fmaxi.

⑫ Adjustment controls **+** and **-**

These two controls allow selected parameters to be increased or decreased after making a selection with the menu key ⑪: (Inspiratory trigger, I/E ratio, Pmini or maxi.) (See Chapter III.3 Adjusting **OSIRIS 2**.)

Changes to settings are only taken into account after the VALID button (14) has been pressed.

(13) Screen highlight button "☀"

This button allows highlighting of the screen (1) for better visibility of information when using the equipment at night or in low light levels. Each time the button is pressed the highlight remains on (for 30 seconds if running on batteries; until the next time the button is pressed if running off the mains). If battery charge is low this facility will not be available.

(14) "→" button

This button validates changed ventilation parameter values.

(15) "DATA" button, for displaying current patient parameters (Frequency and I/E ratio)

II. 2. Electrical and pneumatic connections on side of unit

(16) Oxygen supply

Foolproofed connector according to the national regulations in force: From 280 kPa (2.8 bar) to 600 kPa (6 bar).

(17) Patient circuit connection

Output cone standardized to ISO 22/15.


(18) Expiratory valve control connection

This is the connection for the tube that pilots inflation and deflation of the expiratory valve diaphragm. It controls the PEEP level during the expiratory phase.

(19) External electrical supply connection

Made with a JACK type connector.

When used with a patient who is in recovery room, **Air Liquide Medical Systems** strongly recommends using the cable clamp (24) on the side of the equipment to ensure optimum safety against accidental unplugging. The cable clamp must be fixed on the rear part of the unit approximately 5 cm from the electrical outlet.

II. 3. Rear of unit

⑳ Ambient air intake and battery access cover

Incoming air is filtered by a foam pad.

The air intake must remain totally unobstructed.

The housing ⑳ contains a battery pack, NiCd 700 mAh type (8 elements of 1.2 V in series).

OSIRIS 2 can function on this internal battery pack or from the mains supply. The battery pack charges when the equipment is connected to the mains supply ⑲ via the special electrical supply unit.

Use only the mains adapter and the battery pack specified and supplied by Air Liquide Medical Systems. Using any other adapter or battery pack releases Air Liquide Medical Systems free of any responsibility.

A new battery will only provide maximum performance when it has been charged and discharged two or three times!


㉑ Manufacturer's plate

This plate indicates the serial number of the equipment, and CE_{0459} indicates that **OSIRIS 2** conforms to European directive 93/42/CEE for medical apparatus.

㉒ Attachment feet (for rails, beds, etc.).

㉓ Carrying handle

㉔ Label giving instructions for setting-up **OSIRIS 2**

III. SETTING UP THE UNIT

III. 1. Installing the battery

OSIRIS units are supplied with the battery not installed.

To install the battery:

- Remove the cover plate (20),
- Connect the battery and insert it in its housing,
- Reinstall the cover plate, ensuring that it is fully engaged.

III. 2. Connections

Oxygen supply to OSIRIS 2

- Connect the **OSIRIS 2** input to a **source of oxygen** via the O₂ supply tube.

If the oxygen is from a cylinder, the cylinder must be equipped with a reducer adapted to supplying medical oxygen.

The nominal supply of O₂ must be between 280 kPa (2.8 bar) and 600 kPa (6 bar).

⚠ Connect the gas supply hose to the ventilator connector before connecting it on the network (or cylinder) side. To disconnect (after shutting down the unit), start by disconnecting the hose from the network (or cylinder) connector.

Electrical supply connection

- When running off the mains supply, connect the special electrical unit (supplied by **Air Liquide Medical Systems**).

The electrical powersupply also charges the internal battery, whether OSIRIS 2 is in use or not (battery charging time: 10 hours).

Setting up the patient circuit

Two types of patient circuit can be used on the **OSIRIS** range.

Single-use patient circuit


A single-use patient circuit, composed of:

- a single-use patient tube (1) together with its control tube (2),
- a single-use expiratory valve already connected to the tubes above (3).

All our single use patient circuits are assembled in their original packaging.

Setting up

- Connect one end of the ringed tube onto the outlet of the **OSIRIS 2** (standard conical outlet, 22 mm diameter).
- **Whenever possible, fit a new bacteriological filter onto the "patient" output of the expiratory valve.**

The expiratory valve (ref. KY133900) is for use with OSIRIS 2 exclusively.

Autoclavable patient circuit (optional)

An autoclavable patient circuit, composed of:

- A reusable patient tube (1) (blue silicone) together with its control tube (2),
- A reusable expiratory valve (3).

To carry out maintenance procedures on the patient circuit, refer to the user manual (§IV.4 for OSIRIS 2)

Setting up

- Connect one end of the ringed tube onto the outlet of the OSIRIS 2 (standard conical outlet, 22 mm diameter).
- Whenever possible, fit a new bacteriological filter onto the "patient" output of the expiratory valve.
- Connect the other end of the ringed tube to the input of the expiratory valve (see schema below).
- Connect the OSIRIS 2 control to the diaphragm of the expiratory valve via the small diameter open-ended tube. When the connections are finished check that the small diameter tube is not pinched.

Starting up

- Switch to the operating mode required; VAC/ACMV 100 % or VAC/ACMV Air + O₂ or VS-AI/PS.

Check that the alarm circuit is working when starting up OSIRIS 2. All the lamps should flash and an audible signal should sound.

- Set the Pmaxi alarm threshold. Then close off the patient outlet and check visually on the bargraph that the patient circuit pressure is limited to (Pmax) (Pmax adjustable from 10 to 80 cm H₂O).
- In VS-AI/PS mode, check that the selected pressure value (2nd line of the Pmax display) is close to the peak pressure displayed on the bargraph to be sure that the sensor is working correctly.
- Check that everything is working normally (see checklist before use).

The unit is ready and may be connected to the patient.

- Fix a mask, or intubation probe, or tracheostomy cannula to the expiratory valve (or to the « \updownarrow » connection of the expiratory valve if you are using an autoclavable patient circuit (diagram opposite))
- **Whenever possible, fit a new bacteriological filter onto the "patient" output of the expiratory valve.**

If the case of any abnormal functioning of OSIRIS 2, refer to the paragraph "Operating incidents"

III. 3. Adjusting OSIRIS 2

Choice of ventilation mode

The mode selector allows **OSIRIS 1 en** to be started up and for the ventilation mode most suited to the patient's needs to be selected: VC-VAC 100% O₂ mode, VC-VAC Air + O₂ mode, VS-AI mode.

When changing from the VS-AI/PS mode to VAC/ACMV Air + O₂ mode or VAC/ACMV 100% O₂ mode, the VT setting may be critical. A screen message "Adjust VT ?" will appear, plus an alarm. If an adjustment is not made within 7 seconds insufflation restarts with the displayed value of the previous setting.

Mode Volume Assisté Contrôlé (VC/VAC) Assisted Controlled Ventilation mode (VAC/ACMV)

The VAC/ACMV 100% and VAC/ACMV Air + O₂ modes are controlled and assisted controlled modes.

- In **VAC/ACMV 100% O₂** mode the FiO₂ supplied is 100 %.
- In **VAC/ACMV 100% Air + O₂** mode the FiO₂ supplied is about 60%.

A standard tidal volume is fed into the patient circuit on each inspiratory effort of the patient helped by a constant flowrate.

If there is no patient activity, the equipment ventilates automatically on set parameters (F mini, etc.).

Because the unit works with a "venturi" system, the FiO₂ supplied in the Air + O₂ position depends on patient ventilation:

- if $\dot{V} = 5$ L/min : 80 % +/- 10 % of O₂
- if $\dot{V} = 10$ L/min : 70 % +/- 10 % of O₂
- if $\dot{V} = 20$ L/min : 60 % +/- 10 % of O₂
- if $\dot{V} = 30$ L/min : 60 % +/- 10 % of O₂

As soon as the patient triggers a respiratory cycle, the « SPONT CYCLE » message is briefly displayed.

Setting the frequency

Frequency is set by turning the graduated knob ⑨.

Frequency can be set from 6 to 40 cycles/min.

This setting directly affects the V_T insufflated to the patient. It is therefore necessary to reset the V_T after each respiratory frequency setting, using the inspiratory flowrate setting knob ⑩.

Setting V_T ((VAC/ACMV 100% mode and VAC/ACMV Air + O₂ mode)



This setting is made by turning the inspiratory flowrate knob ⑩.


The setting range for V_T is a function of the inspiratory time setting (that is, the frequency and I/E ratio chosen).

V_T is continuously adjustable from 100 ml to more than 3000 ml.

The value of the VT parameter is displayed in the lower part of the **OSIRIS 2** screen.

When we switch from 100 % O₂ VAC mode to Air O₂ VAC, the current volume of the instruction displayed on the screen is slightly modified.

Therefore, the position of the settings button must be adjusted when we switch from one mode to another.

Setting the PEEP (Positive End Expiratory Pressure)

This setting is made by turning the graduated knob (3). The PEEP is adjustable from 0 to 15 cmH₂O.

A mechanical safety device (25) has been fitted to this knob to avoid any possible setting of PEEP > 10 cmH₂O. To set a PEEP higher than this value, press on the metal finger next to the knob.

Setting other ventilation parameters (inspiratory trigger, I/E ratio, Pmini alarm, Fmaxi alarm)

Press successively on the (3) button to display the name of the parameter to be set, on the first line of the screen:

- Inspiratory threshold trigger (TRIG)
 - Can be blocked (value = OFF)
 - Keeps its value between two « power on » switches
- I/E ratio (I/E)
- Pmini alarm threshold (Pmini)
- High frequency alarm threshold (Fmaxi).

Then set the parameter to the desired value using the (+) and (-) buttons.

Validate by pressing on the (↻) button and the setting will be stored by **OSIRIS 2**.

Repeat the operation for each parameter.

The I/E ratio setting directly influences the V_T insufflated to the patient. Therefore, after changing the I/E ratio, it may be necessary to reset the V_T using the inspiratory flowrate setting knob (10).

Setting the Pmaxi safety value

This setting is made by turning the graduated knob (4). Pmaxi can be set between 10 and 80 cmH₂O/hPa. The setting is indicated on the screen by a vertical line.

Each time the device is switched on the default values for the settings parameters are used: TRIG = -1 cmH₂O, I/E = 1/2, RRmax = 40 bpm, Pmin = 5 cmH₂O for software version V1.020 or earlier and 10 cmH₂O for software version V1.021 or later. To find out the software version installed, switch on **OSIRIS 2**. The message «OSIRIS VX.XXX» will appear during startup. The version number is X.XXX.

Pressure Support Ventilation (PSV)

PSV is a **pressure-assisted** spontaneous ventilation mode.

A constant positive pressure is maintained in the patient circuit with each inspiratory effort by the patient. The support flow is generated in order to maintain constant inspiratory pressure.

When the pressure support setting is less than the PEEP, a message appears on the screen. The pressure support must then be increased.

The inspiratory phase ends when an expiratory effort from the patient is detected, or the TI setting is reached.

FiO₂ in PSV mode depends on the ventilation parameters, resistance, and patient compliance. It remains between 65 and 95%.

 *When switching from PSV mode to VCV/ACMV Air + O2 or VCV/ACMV 100%, the VT setting can be important. A message 'Adjust VT?' appears on screen with an alarm. If the adjustment is not made within seven seconds, insufflation is resumed with the displayed VT value, which is based on the previous setting.*

 *As soon as the patient triggers a respiratory cycle, the message "CYCLE SPONT" is briefly displayed.*

 *When the respiratory cycle is initiated by the machine, the message "CYCLE CONTROL" is briefly displayed.*

Setting the Minimum Respiratory Rate

This is a continuous adjustment, performed by rotating the graduated knob **9**.

The frequency can be adjusted from 6 to 40 bpm.

 *In case of patient inactivity, the device automatically applies pressure-controlled ventilation according to the set parameters (F min, I/E, etc.).*

 *As stated in the 'expiratory trigger' paragraph, the Fmin setting affects the duration of the refractory period of the expiratory trigger.*

Adjusting the PEEP (Positive End-of-Expiration Pressure)

This adjustment is performed by rotating the graduated knob **3**.

The PEEP setting is adjusted continuously from 0 to 15 cmH₂O.

Note: A mechanical safety device is placed on this knob to prevent any accidental setting of PEEP > 10 cmH₂O.

To set a higher PEEP value, press the metal index pin located beside the knob.

Adjusting the pressure support

This adjustment is performed by rotating the inspiratory flow rate knob **10**.

The pressure is adjustable from 3 to 44 cmH₂O.

The value appearing on the second line is the pressure setting (PI max = pressure support + PEEP).

 *When the pressure support setting is less than the PEEP, a message appears on the screen. The pressure support must then be increased.*

Adjusting the other ventilation parameters (inspiratory trigger, I/E ratio, Pmin alarm, Fmax alarm)

Press the MENU key to access each of the settings to be adjusted in turn. The value appears on the first line of the screen:

- Inspiratory trigger threshold (TRIG.);
- I/E ratio (I/E);
- Pmin alarm threshold (Pmin);
- High frequency alarm threshold (Fmax).

Next, adjust the parameter displayed to the desired value using the keys **+** or **-**.

Confirm by pressing the key. OSIRIS 2 will then take the setting into account.

Repeat the operation for each parameter.

 *As stated in the 'Expiratory trigger' section, the I/E ratio setting affects the duration of the refractory period of the expiratory trigger.*

Adjusting the Pmax safety setting

This is a continuous adjustment, performed by rotating the graduated knob **4**.

Pmax is adjustable from 10 to 80 cmH₂O (hPa).

The setting is displayed directly on the screen of the device, represented by a vertical line.

Expiratory trigger

The expiratory demand is based on the detection of pressure variations during the inspiratory phase.

At the beginning of insufflation, **OSIRIS 2** has a refractory period, during which the expiratory trigger may not be triggered.

The duration of this refractory period varies according to the duration of the inspiratory time: it is equal to half of the inspiratory time setting, with a minimum limit of 300 ms and a maximum limit of 1000 ms.

To limit any incompatibility between patient and machine, the minimum frequency setting should be slightly less than that of the patient (approximately 20% less), which will give the best possible adjustment of the refractory period.

This refractory period, moreover, has the advantage of ensuring at least a certain minimum inspiratory time, and therefore at least a certain minimum tidal volume for any patient who might tend to cut his or her inspiratory time too short.

If the practitioner wishes to place top priority on this safety feature, note that the refractory period may be deliberately extended (although only up to one second) by lowering the frequency setting or by increasing the value of the I/E ratio.

This safety feature is useful in an emergency and transport ventilator, where monitoring is bound to be more limited than on a resuscitation ventilator.

As an indication: on **OSIRIS 2**, this refractory period is at its maximum (1000 ms) for any frequency setting less than or equal to 10 bpm (case I/E = 1/2), and for any frequency setting less than or equal to 15 bpm (case I/E=1/1).

 Each time the device is powered up, the default settings are used: TRIG = -1 cmH₂O, I/E = 1/2, Fmax = 40 bpm, Pmin = 5 cmH₂O for software versions V1.020 or earlier, and 10 cmH₂O for software versions V1.021 or later.

III. 4. Display of frequency and I/E ratio measured by **OSIRIS 2**

OSIRIS 2 displays measured values of the I/E ratio and respiratory frequency **du patient** when the button is pressed successively.  Parameters are displayed in the upper part of the screen for 5 seconds.

 The respiratory frequency measured by **OSIRIS 2** can, in certain cases (respiratory cycles started by the patient), be greater than the frequency set on the unit.

V. DO NOT FORGET

IV. 1. Safety and protection

IV.1.1. Additional air inlet

A valve allows the patient:

- to breath ambient air via **OSIRIS 2** if the unit is stopped (voluntarily or accidentally),
- and, to breath ambient air if the ventilation flowrate is insufficient.

IV.1.2. Pmaxi safety

Pmaxi safety allows pressure limitation and therefore avoids any barotrauma (inspiratory phase stopped, proceeding to expiratory phase when the set threshold is exceeded). This can be set to between 10 and 80 cm H₂O.


IV.1.3. Pmax safety valve

This limits the pressure in the patient circuit to 90 mbar (cmH₂O) in the event of a Pmaxi safety failure.


IV.1.4. Gas filtering

This is assured by:

- a foam pad on the ambient air intake.

- a nylon filter on the O₂ supply.

 A bacteriological filter (type Pall BB50) can be fitted either on the patient circuit output connection or on the "patient" outlet of the expiratory valve.

IV. 2. Alarms

Pmini alarm: Intermittent audible and visual alarms will be activated if pressure in the patient circuit does not reach the set Pmini value (from 5 to 10 mbar (cmH₂O)) during 2 respiratory cycles.

The audible alarm may be inhibited for 2 minutes by pressing on the  button.


Failure/Message	Audible signal	Probable cause	Proposed corrective action
F max ALARM	Intermittent	The respiratory rate measured by OSIRIS is higher than the Fmax threshold set by the user	Increase the value of the alarm threshold or resolve the problem of tachypnea
		Autotriggering of the inspiratory trigger	Decrease the sensitivity of the inspiratory trigger
PI < Pmin ALARM	Intermittent	The pressure in the patient circuit does not exceed the threshold (PEEP + Pmin) set by the user for inspiration, for at least 2 respiratory cycles	Decrease the Pmin threshold, or increase the set VT, or the set PI.
		Patient circuit leak	Ensure the connections are leaktight
		Leak at expiratory valve level	Fit the diaphragm correctly or replace it
		Probe's cuff not leaktight	Inflate the probe's cuff and check that it is leaktight
		No O ₂ supply, O ₂ cylinder empty	Connect the O ₂ or replace the O ₂ cylinder
		Expiratory valve control tube clamped	Check the expiratory valve control tube
Pmax ALARM	BEEP	The pressure in the patient circuit has reached the set Pmax value, for one or two respiratory cycles	
Pmax ALARME	Intermittent	The pressure in the patient circuit has reached the set Pmax threshold value, for at least three respiratory cycles	Increase the alarm threshold value or decrease the set VT or the set PI
Machine anomaly. Ventilation shutdown. Patient exposed to ambient air	Continuous	Electronic failure	Contact the technical department
P. max anomaly	Intermittent	Electronic failure	Ventilation can continue but without pressure monitoring Monitor the patient closely. Contact the technical department

Failure/Message	Audible signal	Probable cause	Proposed corrective action
FLAT BATTERY	Intermittent	Flat battery	Ventilation can continue but the unit may shut down at any time. Switch immediately to the electrical power supply or install a charged battery.
PEEP > 20 cmH ₂ O	Intermittent	Expiratory valve incompatible with the respirator	Ensure that the expiratory valve used is compatible with OSIRIS (ask us)
		Patient valve outlet obstructed	Check that nothing is hindering the patient's expiration
ONE BUZZER OUT OF SERVICE	Intermittent	One of the two buzzers is damaged	Contact the technical department
CHECK SENSOR	Intermittent	1 or 2 pressure sensors damaged	Contact the technical department
DISCONNECTION	Intermittent	The patient has been disconnected	Reconnect the patient
		Incorrect connection of one of the elements of the patient circuit	Ensure that the connections are leaktight
CLOCK BATTERY OUT OF SERVICE	NO	Lithium battery nearly flat	Ventilation continues. Risk of incorrect event dating in black box. Contact the technical department
CALIBRATION PROBLEM	Intermittent	Calibration failure during ventilation	Ventilation continues with the default parameters. Contact the technical department
BUZZERS OUT OF SERVICE	NO	Damaged buzzers	Ventilation continues without audible alarm Contact the technical department
Start-up problem		Flat battery	Connect to the electrical power supply or install a charged battery
		Unit failure	Contact the technical department
Minute ventilation less than the set minute ventilation		Ambient air inlet obstructed	Release the ambient air inlet Replace the foam filter
LED is flashing		Battery missing or out of service	Install a new battery

IV. 3. Operating incidents

OBSERVATIONS	PROBABLE CAUSES	CORRECTIVE ACTION
No ventilation	<ul style="list-style-type: none"> • Patient disconnected * • Leak in the patient circuit * • Control pipe clamped * • No oxygen supply * 	<ul style="list-style-type: none"> • Check the seal of the mask or probe • Check the patient circuit • Check the control tube • (its fitting, wear and tear, etc.) • Check the connections and the oxygen source
Start-up problem	<ul style="list-style-type: none"> • Mode selector in Off position • Microprocessor board battery out of service • Device fault 	<ul style="list-style-type: none"> • Switch the mode selector to one of the 3 operating positions • Contact A.S.D • Contact A.S.D
No start-up with audible and visual signal	<ul style="list-style-type: none"> • Battery flat 	<ul style="list-style-type: none"> • Change the battery or recharge it.
Ventilation still insufficient	<ul style="list-style-type: none"> • Ambient air inlet obstructed • Leak in patient circuit * 	<ul style="list-style-type: none"> • Clear the air inlet • Replace the foam filter • Check the patient circuit
Pmin <ul style="list-style-type: none"> • The battery LED is flashing • The Pmin alarm is permanently triggered 	<ul style="list-style-type: none"> • Battery failure if the device is used with internal batteries • Pmin incorrectly set 	<ul style="list-style-type: none"> • Charge the batteries, change them if they are defective • Reset the Pmin threshold

* The fault is signaled by the PMINI alarm

Material		Autoclave 136°C 220kPa 20 minutes	Ethylene oxide 55°C
Silicone 		OK	OK
Polysulphone 		OK	OK
PVC 		No	OK

IV. 4. Cleaning, Disinfection, Sterilization

IV.4.1. Single-use patient circuit

A single-use patient circuit, which is composed of:

- A single-use patient tube (1) together with its control tube (2);
- A single-use expiratory valve already connected to the tubes above (3).

 All our single use patient circuits are assembled in their original packaging.

 Three single-use patient circuits (KG019900) are delivered together with the ventilators.

IV.4.2. Autoclavable patient circuit (optional)

- Clean the patient circuit with soapy water or with a cleaning solution..

 Do not use powerful solvents (ether, acetone, etc.).

- Disinfection of the patient circuit may be made by immersion in a disinfecting solution, **after cleaning**.
- Steam sterilize the patient circuit, following the instructions in the adjacent table.
- Bacteriological filter: It is advisable to change the bacteriological filter on the insufflation output connection at least each time the patient circuit is sterilized, or each time the equipment is connected to a new patient.

IV.4.3. OSIRIS 2

- The unit may be cleaned with a cloth dampened with soapy water, avoiding any risk of water entering the apparatus (electronic circuitry may be damaged), then wipe with a dry cloth or a cloth impregnated with a water/ alcohol solution.
- Disinfection of the outer surfaces of the unit may be made with a spray such as Surfanios® (adhering closely to the manufacturer's instructions, especially dosage) or with cloths impregnated with an alcohol-based disinfecting solution.
- **OSIRIS 2** may be sterilized with a formaldehyde gas unit.

Procede as follows:

- Connect the **OSIRIS 2** to the sterilizing unit via the ambient air connection (20) on the rear of the ventilator,
- Set button f to 20 c/min,
- Run the sterilizing unit according to its manufacturer's instructions.

It is imperative to use air which is extremely clean and free of any trace of oil or water.

Sterilization with formol using any other procedure than that described above may cause formol crystallization and stop OSIRIS 2 working.


IV. 5. Maintenance

IV.5.1. By the user

- Before each use, check the patient circuit (**expiratory valve, expiratory diaphragm valve, etc.**), the oxygen supply tube, and the general condition of the external electrical cable...

Also check that the "X" indicator is not lit.

- Replacing the battery-pack
 1. Remove the cover;
 2. Remove the battery pack and its connector,

3. Disconnect the battery pack and replace it with a new (or charged) battery pack,
4. Replace the battery pack and its connector in its housing,
5. Replace the cover, pressing until it clicks shut.

 If the unit is not running off the mains supply when the battery pack is changed, **OSIRIS 2** restarts with these default values: $TRIG = -1 \text{ cmH}_2\text{O}$, $I/E = 1/2$, $RR_{max} = 40 \text{ bpm}$, $P_{min} = 5 \text{ cmH}_2\text{O}$ for software version V1.020 or earlier and $10 \text{ cmH}_2\text{O}$ for software version V1.021 or later.

IV.5.2. Prolonged storage or unit not in use

If your **OSIRIS** ventilator will remain unused for a period exceeding one month, its battery should be removed and placed in storage according to the following recommendations:

- Charge the battery fully before storage.
- Recharge the battery at least once every six months.
- Store the battery in a location at moderate temperature. Heat tends to increase the rate of discharge, as shown in the table below:

Storage temperature	Battery self-discharge rate
0 °C	10 % / month
20 °C	20 % / month
40 °C	40 % / month

If completely discharged, the battery can be reactivated by applying the following procedure:

- If the battery lost its charge whilst connected to your **OSIRIS** ventilator, run four or five complete charge/discharge cycles.
- If the battery lost its charge during storage (disconnected), perform a full charge.
- In all cases, check battery life after reactivation under the following ventilation conditions: ACMV 100%, frequency = 20 c/min, I/E ratio = 1/2 and tidal volume = 500 ml. You should obtain 6 to 12 hours of battery operation.



 When the battery is completely discharged, you can always operate your **OSIRIS** ventilator using the charging cable (supplied).

IV.5.3. By technicians

These operations must be made by technicians who have received the appropriate training.

Every year

- Replace:
 - The various filters of the unit (foam filter on ambient air intake ②, O₂ filter ①, etc.)
 - The internal battery pack,

- Carry out a complete check on **OSIRIS 2**

For more details refer to the Maintenance manual.

IV. 6. Accessories

The accessories used with **OSIRIS 2** must:

- be oxygen compatible,
- be biocompatibles,
- conform to the general requirements of EN 60601-1 and of directive 93/42/CEE,
- and not be antistatic or conductors of electricity.

The use of accessories which do not conform to the above descriptions releases **Air Liquide Medical Systems** from any responsibility in the event of an incident.

Accessories referenced by **Air Liquide Medical Systems** or included in the batch of accessories delivered with the equipment conform to these requirements.

Electromedical apparatus associated or used with **OSIRIS 2** must conform to the requirements of the directive 93/42/CEE (medical equipment).

The apparatus is delivered with:

1 transportation case	KF006500
3 single-use patient circuits	
1 NiMH battery pack	KY564600
1 mains/charger adapter	KY564800
1 adult mask	KM030000
1 child's mask	KM030100
1 3 m O ₂ tube with fixed connector	BF030200

Optional accessories:

1 stand for OSIRIS 2	KA004600
1 autoclavable patient circuit	KG018700
1 autoclavable expiratory valve	KY133900
20 single-use patient circuits	KG020300
1 articulated arm, chrome	KB005200
1 wheeled base	KB005000
1 Ambulance compatible wall-support	KY581900

IV. 7. Waste disposal

All waste products coming from **OSIRIS 2** (the patient circuit, etc.) must be disposed of according to the appropriate procedures of the hospital.

IV. 8. Disposal of equipment

To protect the environment, all equipment disposal must be made according to the appropriate procedures of the hospital. Besides, traceability being imposed by the **CE** marking, it is obligatory to communicate the disposed equipment's serial number to **Air Liquide Medical Systems** technical department.

V. TECHNICAL DESCRIPTION

V. 1. Operational principles

OSIRIS 2 is a frequency controlled ventilator (from 6 to 40 cycles/min), with adjustable I/E ratio (from 1/1 to 1/3 by steps of 0.5.).

The gas supplied may be pure oxygen or a mixture of ambient air/oxygen (FiO_2 at about 60%).

The PEEP can be adjusted from 0 to 15 hPa (cmH_2O). An adjustable Inspiratory Trigger level from -4 to -0.5 cmH_2O (in steps of 0.5 cmH_2O) assures maximum patient comfort.

An adjustable minimum pressure alarm (5 to 10 cmH_2O above the PEEP level, in steps of 1 cmH_2O), as well as an adjustable Fmaxi alarm (from 20 to 60 cpm, in steps of 5 cpm) and an adjustable Pmaxi safety alarm (from 10 to 80 cmH_2O) guaranteeing sure and efficient ventilation

The device is also fitted with:

- an ambient air inlet valve which allows the patient to breathe freely through the device should a breakdown occur,
- a safety valve calibrated to 90 mbar (cmH_2O).

V. 2. Technical characteristics

- Pressure displayed by bargraph of - 4 to 80 cmH_2O ($\pm 1 \text{ cmH}_2\text{O}$).
- Pneumatic supply: from 280 kPa (2.8 bar) to 600 kPa (6 bar).
- Electrical supply for alarms: battery pack, type NiMH 700 mAh. (8 elements of 1.2 Volt in series).
- Battery autonomy: 12 hours (VAC/ACMV 100 %, frequency = 20 c/min, I/E ratio = 1/2 and $V_T = 500 \text{ ml}$).
- Input: 230 V ~ / 50 Hz Output 15 V --- 6 W or 110 V ~ / 60 Hz Output 15 V --- 6 W (consult us).
- Electrical class: Class II
- Type of apparatus: Type B
- Power loss protection: Assured by the "Mains supply fault" alarm.
- Consumption in use:

- The unit alone: 1 l/min.
- In total: patient consumption + consumption of unit.

Example of consumption:

- For an adult, the average consumption is 6 l/min of gas (air or oxygen).
- The consumption of the unit is fixed at 1 l/min..
- In the case of the use of a type B compact cylinder, the cylinder volume can be 5 litres. The gas being compressed to 200 bar, there is, in this case, 1000 litres of gas available.
- In the case of our example we have available therefore an autonomy of about 2 hours and 20 minutes if the ventilation is made at a concentration of 100 % oxygen.
- Inspiratory and expiratory resistance of complete system: (reusable circuit without bacteriological filter):

Flowrate	Insp R (cmH_2O)	Exp R (cmH_2O)
30 L/min	1.6	0.2
60 L/min	3.4	0.5






- Flowrate and volume measurements are given in ATPD.
- Mode of operation: Continuous service
- Operating temperature limits: from -18°C to $+50^\circ\text{C}$ (relative humidity < 95 %).
- Storage temperature: from -40°C to $+60^\circ\text{C}$ (relative humidity < 95 %).
- For use in atmospheric pressure: from 700 to 1060 mbar.
- Size L x W x H = 250 x 210 x 170 mm
- Weight: 5.0 kg
- Volume of patient circuit with a standard 1500 mm tube: 570 ml for the tube + 10 ml for the valve.

V. 3. REGULATIONS AND STANDARDS

- Directive 93/42 EEC concerning medical devices
- Directive 95/54/CE
- Standard EN 60 601-1 and its collateral standards concerning the safety of electromedical devices.
- Standard EN 794-3 and its collateral standards concerning transport and emergency ventilators.
- EN 1789: Medical transport vehicles and their equipment

Year of **CE**-marking : 1999

V. 4. Symbols and abbreviations

$P_{EHP/PEEP}$ (cmH ₂ O/hPa)	Positive End Expiratory Pressure.		Setting button: increases the selected parameter value.
A.S.S.	After Sales Service.		Setting button: reduces the selected parameter value.
I/E	Inspiratory time/Expiratory time ratio		Validation button: validates a changed ventilation parameter value.
F/RR (cpm/bpm)	Respiratory frequency (cycles per minute).		Button for displaying patient parameters during ventilation (f and I/E).
P (cmH ₂ O/hPa)	Pressure value (in cmH ₂ O or hPa)		Screen highlighting
	Important: Consult accompanying documents.		Audible alarm inhibit
	Type B apparatus	ALARM	Visual alarm indicator
	Class II apparatus (double insulation)		Mains supply operation indicator
	Conforms to directive 93/42/CEE (established by organization No. 0459).		Battery supply operation indicator
IP34	Protected against penetration by solids of diameter greater than or equal to 2.5 mm and against water-spray ingress.		Limited battery autonomy indicator (≤ 60 minutes)
	Do not grease.		Oxygen supply (③ : fixed O ₂ outlet)
	Fire risk. Keep aware from all inflammable materials or liquids.		Patient connection (① : patient circuit; ② : expiratory control valve circuit)
P_{maxi} (cmH ₂ O/hPa)	Maximum operating pressure setting (in cmH ₂ O or hPa).		Electrical supply
0	Stop		Sense of gas circulation
$V_{AC}/ACMV$ 100% O ₂	Ventilation in VC - VAC mode with an FiO ₂ of 100 %.		see § concerning internal electrical supply.
$V_{AC}/ACMV$ A ₁ + O ₂	Ventilation in VC - VAC mode with an FiO ₂ > 50 % (O ₂ / ambient air mixture).	e2	Compliance with Directive 72/245/EEC and 95/54/EC under homologation number 03090.
	Selection button for parameters to set (other than V _T , frequency and PEEP), that is: I/E, inspiratory trigger, P _{mini} and F _{maxi} .		

Material in direct or indirect contact with the patient

Reducer:

- Technopolymer
- Nitrile
- Brass
- Phosphated steel
- Stainless steel

Flowrate valve:

- Steel
- Nitrile
- Stainless steel
- Aluminium
- Brass
- Delerin
- EPDM

Safety valve:

- Aluminium
- Silicone
- Stainless steel
- Brass

Electro-distributor:

- Aluminium
- Stainless steel
- POM
- Nitrile

Output connector assembly:

- Nitrile
- Stainless steel
- Silicone

Air intake assembly:


- ABS - VO
- Bulpren
- PVC

Tubing:

- Rilsan
- Medical polysulphone
- Polyurethane

 **OSIRIS 2** and its associated accessories do not contain any Latex.

Pneumatic operation schematic

The details given in this section are for information purposes only.

In the interest of constantly improving our products, we reserve the right to incorporate any modification considered useful.

CHECK-LIST: Before use / Return from maintenance

	Yes Done	No Not done
<ul style="list-style-type: none"> • Check that the mains supply cable and the pipes (patient circuit, expiratory valve, expiratory valve control tube and the oxygen supply pipe, etc.) are in a satisfactory general condition. 	<input type="checkbox"/>	<input type="checkbox"/>
<ul style="list-style-type: none"> • Connect a gas supply to OSIRIS 2 from an oxygen source via an oxygen supply tube. Then check that the supply pressure is between 280 kPa (2.8 bar) and 600 kPa (6 bar) inclusively using a pressure monitoring system or an independent control pressure gauge. 	<input type="checkbox"/>	<input type="checkbox"/>
<ul style="list-style-type: none"> • If using an oxygen cylinder, this must be equipped with reducer. Check that the supply pressure is between 280 kPa (2.8 bar) and 600 kPa (6 bar). 	<input type="checkbox"/>	<input type="checkbox"/>
<ul style="list-style-type: none"> • Connect the patient circuit (cf Pages 7 and 8). 	<input type="checkbox"/>	<input type="checkbox"/>
<ul style="list-style-type: none"> • Switch the device on (switch in VAC/ACMV 100%, VAC/ACMV Air + O2 or VS-AI/PS). An audible and visible signal will be activated. Check that the "⊗" LED is not lit. 	<input type="checkbox"/>	<input type="checkbox"/>
<ul style="list-style-type: none"> • Open the patient output and check that the pressure read on the bargraph is lower than 10 cmH₂O. The Pmini alarm must activate after 2 cycles have occurred. 	<input type="checkbox"/>	<input type="checkbox"/>
<ul style="list-style-type: none"> • Close off the patient outlet, and check that the PMAX audible and visual alarm is activated. 	<input type="checkbox"/>	<input type="checkbox"/>
<p>Disconnect the electrical power-supply. Check that the unit switches over to battery supply without stopping ventilation.</p>	<input type="checkbox"/>	<input type="checkbox"/>
<ul style="list-style-type: none"> • Connect OSIRIS 2 to the patient. The apparatus is ready to be used. 	<input type="checkbox"/>	<input type="checkbox"/>
	<input type="checkbox"/>	<input type="checkbox"/>
	<input type="checkbox"/>	<input type="checkbox"/>

1st year

Date:

Number of hours:

Technician's name:

.....

Signature
and
stamp:

2nd year

Date:

Number of hours:

Technician's name:

.....

Signature
and
stamp:

3rd year

Complete maintenance

Date:

Number of hours:

Technician's name:

.....

Signature
and
stamp:

4th year

Date:

Number of hours:

Technician's name:

.....

Signature
and
stamp:

5th year

Date:

Number of hours:

Technician's name:

.....

Signature
and
stamp:

6th year

Complete ADR maintenance

Date:

Number of hours:

Technician's name:

.....

Signature
and
stamp:

7th year

Date:

Number of hours:

Technician's name:

.....

Signature
and
stamp:

8th year

Date:

Number of hours:

Technician's name:

.....

Signature
and
stamp:

Maintenance form

OSIRIS 2 NO.:

.....

In service on :

.....

.....

Maintenance assured by:

.....

.....

.....

Your distributor:

.....

Address:

.....

.....

.....

.....

.....

.....

.....

Telephone :

.....

Preventive maintenance of equipment must be done according to instructions given by the manufacturer in maintenance manuals and possible updates.

Operations must be performed by technicians who have received appropriate training.

Use only original spart parts.

On request, the supplier will make available: circuit diagrams, component lists, technical descriptions, or any other information useful to qualified technical personnel, to repair parts of the equipment stated as repairable by the manufacturer.

Manufacturer:
Air Liquide Medical Systems
 Parc de Haute Technologie
 6 rue Georges Besse CE 80
 92182 Antony CEDEX - FRANCE
 Tél : 33 (0)1 40 96 66 00
 Fax : 33 (0)1 40 96 67 00

Website:
www.airliquidemedicalsystems.com
 Hotline:
ALmedicalsystems.services@airliquide.com

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Air Liquide Medical Systems S.A.